

Product datasheet for **MR219430**

Fam47e (NM_001033478) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Fam47e (NM_001033478) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Fam47e
Synonyms: C730043O17; Gm1381
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >MR219430 representing NM_001033478
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGCTGTCCTGCATCCAGCCGTTCCACCTTACATTTCTTTGCATCATCCGTGAGCACCTGGCATGTG
ATTGTCTGCCTCGATGTCTCTCCATCCTGCATAGCCGAAGCCACGCAAAGGCCAAGAAGACGCTGAT
CTCACGTTCAAGCCAGTGGCTTTGTGAAGAAAAGTCAAGCAAAGTGGATCACTCTATAAGACAGACTC
CTCCATGATGACGTACACAGAGGAGTACCGACTTCTGCCACTGGCCGAGGACCTTGAAGCTCTGCCA
TTGAGGAAGATTTCGTCCTGCAGCAGTTTACATCGGCTATCAGACTAGACGCAGCTGTGATGCGCTTCT
CAGGCTGAGGCTGAACCAGGACAACTTTCATCAATAGACCACCTTCCCTGCGAGAACACTGTGGGCAC
TGGAGTGAACCTCGTCCCTCTGGAAGAGCAAATCCTCAGAAGCCGAAGCGGGTAAAGATGAGGTACGGAG
CCTGGTACCTGAACACCAGCTTGTGGAAGAGGCAAAGAGCGGATGAACCTCTGGTCGACCCCATGGTCTC
ACACAAAGCTCAAGATCAACCTTTAAAGAGCAGCTTCAGGAGCAGGGGGAGCTGCTTGCAGGCCTGCGA
GGGACAGCAGCCTTTAAGGACTTCACTTAAGCAGGGGCTACAGGATGCCAGGTTCTTGAAGAATAT
ATGCGGAGGAAAAATAAATCTGAAAAATTAAGGCTCCTAAAAAACTACGCAAACAGAGAAGAATCC
CGGGAGCAGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >MR219430 representing NM_001033478
 Red=Cloning site Green=Tags(s)

MLSCIQPLPTLHFLCIIREHLACDCLPRCLSPSCIAEATQRPKKTLSRSGQWLCEEKSSKVDLSLYKDRL
 LHDDVHRGVTDFCHWAEDLGSSAIEEEFVLQQFDIGYQTRRSCDALLRLRLNQDKLFINRPPSLREHCGH
 WSEPRSSGRANPQKPKRVKMRYGAWYLNLSLWKRQRADEPLVDPMVSHKAQDSTFKEQLQEQLAGLR
 GTAAFKDFILSRGYRMPRFLEKIYAEENKSENKAPKKTQTTEKNPGR

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mm9071_g06.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF

ACCN: NM_001033478

ORF Size: 780 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001033478.2](#), [NP_001028650.2](#)

RefSeq Size: 1505 bp


RefSeq ORF: 783 bp

Locus ID: 384198

Cytogenetics: 5 E2

MW: 30.1 kDa

Product images:



Circular map for MR219430