

## Product datasheet for **MR219449L3V**

### Abat (NM\_001170978) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Abat (NM_001170978) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Abat
Synonyms:	9630038C02Rik; AI255750; GABA; Gabaat; Gabat; Gm9851; I54; Laibat; X61497
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001170978
ORF Size:	1332 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR219449).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001170978.1</a> , <a href="#">NP_001164449.1</a>
RefSeq Size:	4485 bp
RefSeq ORF:	1335 bp
Locus ID:	268860
UniProt ID:	<a href="#">P61922</a>
Cytogenetics:	16 A1



[View online »](#)

**Gene Summary:**

The encoded gene product is responsible for catabolism of gamma-aminobutyric acid (GABA), a mostly inhibitory neurotransmitter in the central nervous system, into succinic semialdehyde. Deficiency of this encoded protein includes psychomotor retardation, hypotonia, hyperreflexia, lethargy, refractory seizures, and EEG abnormalities. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jan 2010]