

## Product datasheet for **MR219457L3V**

### **Kcnk13 (NM\_001164426) Mouse Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Kcnk13 (NM_001164426) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Kcnk13
Synonyms:	BB085247; F730021E22Rik; Gm1570; Gm1685
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001164426
ORF Size:	1215 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR219457).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001164426.1</a> , <a href="#">NP_001157898.1</a>
RefSeq Size:	3075 bp
RefSeq ORF:	1218 bp
Locus ID:	217826
UniProt ID:	<a href="#">Q8R1P5</a>
Cytogenetics:	12 E
Gene Summary:	Potassium channel displaying weak inward rectification in symmetrical K(+) solution. [UniProtKB/Swiss-Prot Function]



[View online »](#)