

Product datasheet for **MR219495**

Sik1 (NM_010831) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Sik1 (NM_010831) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Sik1
Synonyms:	Hrt-20; Msk; Sik; Sik-1; Snf1lk
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>MR219495 representing NM_010831
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGTGATCATGTGCGAGTTCAGTGCAGTCCCTCAGGCACTGGCCAGGGCCAGCAGAAGCCCTCCGGG
 TGGGCTTTTACGACGTGGAACGGACCCTGGGCAAAGGCAATTTTGCAGTGGTTAAGCTGGCGCGCACCG
 AGTCACCAAAAACGCAGGTTGCAATAAAAAATAATTGACAAGACACGTTAGATTCTAGCAATCTGGAGAAG
 ATCTACCGGGAGGTCCAGCTCATGAAACTTTTGAACCACCCAAATATCATCAAGCTTTATCAGGTTATGG
 AGACAAAGGATATGCTCTACATTGTACAGAAATTTGCAAAAAATGGAGAAATGTTTGATTATCTGACTTC
 CAACGGGCACTTGAGTGAACGAGGCTAGGCAGAAATTTGGCAGATTCTGTCAGCCGTGGAGTACTGC
 CACAACCACACATTGTCCACCGGACCTCAAGACGGAGAACCTGCTCCTGGACAGCAACATGGATATCA
 AGCTAGCAGATTTTGGATTTGGGAATTTCTACAAGCCAGGCGAGCCTCTGTCCACGTGGTGTGGAGCCC
 CCCGATGCAGCCCTGAAGTCTTCGAGGGGAAGGAGTACGAGGGTCCCAGCTGGACGCTCTGGAGCCTC
 GGTGTGGTGTACGTCTGTGCTGTGGCTCCCTCCCTTCGACGGGCCAACCTGCCTACGCTGAGAC
 AGCGGGTACTGGAGGGCCGCTTCCGCATCCCTTCTTATGTCTCAAGACTGTGAGACGCTGATCAGACG
 CATGCTGGTCTGGACCCCGTAAGCGCATCACCATAGCCAGATCCGCCAGCACCGGTGGATGCAGGCC
 GACCCACTCTCCTGCAGCAGGATGACCCTGCCTTCGACATGCAAGGCTACACCTCCAACCTGGCGGACT
 ACAACGAACAGGTGCTAGGGATCATGCAGGCCCTCGGCATCGACCGGCAGCGGACTATAGAGTCTCTGCA
 GAACAGCAGCTACAACACTTTGCCGCCATTTACTACCTCTACTCGAGCGCCTCAAGGAGCATCGAAGC
 GCCCAGCCCTCATCCCGGCCACCCCTGCACCCACCAGACAGCCAGCTCCGAAGCTCAGACCTCAGCA
 GTCTGGAGGTTCCCAAGAAATTTCCCGTGTGACCCTTCCGGCCCTCTCTGCTGTGCCACAGCCCCA
 GGCTTGGCTCAGTCTGCTCCTGCAGGCTGAGATAGACTGTGATCTCCACAGCTCACTTCAGCCCTTATTA
 TCCCCCTGGATACCAACTGCAGTGGAGTGTCCGGCACCGATCCATCTCCCCCAGCAGTCTGCTGGACA
 CAGCTATCAGCGAGGAGGCCAGGAGGTTCCAGCCTAGAGGAGGAACAGGAGGTCAGGAACCCCTGCC
 CGGAAGCACAGGCCGGAGGCACACATTGGCTGAAGTCTCCACCCATTTCTCCCGCTCAACCCTCCTTGC
 ATAATTGTCTCCTCTGCCACGGCGAGTCCCTCGGAAGGGACCAGCTCCGACAGCTGCCTCCCTTCT
 CTGCAAGTGAAGTCTGCAGGGCTTGGCAGTGGCTGGCCACCCAGGGCTGCTGGGCACCAGCTCTCC
 AGTCAGATTGGCCTCGCCCTTCTGGGATCACAGTCAGCCACCCCTGTGCTCCAGACTCAGGAGGTTCTG
 GGCACAGCCGTTTACCTCCCGTCACTTCCAGGAGGACGGAGAGCGTCTGATACCTCTCTCACTCAGG
 GGCTGAAGGCCTTCCGGCAGCAGCTGAGGAAAAATGCGAGGACCAAGGGTTCTGGGACTGAACAAGAT
 CAAAGGGTTGGCTCGCCAGGTGTGCCAGTCTCCGTCCGAACCTCCCGGGGAGGGATGAGTACCTCCAC
 ACCCCAGCCCAAGCTCCGGCTGCAAGGCTGCACGACCAGCAATCGGGAGGGCAGGAGCCTACTTGAAG
 AGGTCCTGCACCAGCAGAGGCTGCTCCAGTTACAGCACCCTTCTTACC GCCCGCCCTCATCTGGCTG
 CCAGCAGGGCCCCAGCTGTCCCGGTCCCCTATGTGCTTGTCCCTGTGACAGCCTCCTCGTGTCTGGA
 ATTCCACTGTGCCAACCCCTTCTCCAGGCTGGCATGTCCCTGTGGCGTCTGCTGCGCATCTCCTGG
 ATACCCACTTGACATCAGTGTGGCCAGTGGCCCTCCCCTGAGGGCCCTGCCACAGTGCCTCACCAG
 GCTGTCCCAGGCTGTGATCCAGCTGGGCTGCCACAGGGGACTGTGAGATGGAGGACCTGACCTCTGGC
 CAGCGGGGACATTCGTCTGGTACAG

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA


Protein Sequence: >MR219495 representing NM_010831
Red=Cloning site Green=Tags(s)

MVIMSEFSAVPSGTGGQQKPLRVGFYDVERTLGKGNFAVVKLARHRVTKTQVAIKIIDKTRLDSSNLEK
IYREVQLMKLLNHPNIIKLYQVMETKDMLYIVTEFAKNGEMFDYLTSNGHLSENEARQKFWQILSAVEYC
HNHHIVHRDLKTENLLSDNMDIKLADFGFGNFYKPGPEPLSTWCGSPPYAAPEVFEGKEYEGPQLDVWSL
GVVLYVLVCGSLPFDGPNLPTLRQRVLEGRFRIPFFMSQDCETLIRRMLVVDPAKRITIAQIRQHRWMA
DPTLLQQDDPAFDMQGYTSNLGDYNEQVLGIMQALGIDRQRTIESLQNSSYNHFAAIYYLLERLKEHRS
AQPSSRPTPAPTRQPQLRSSDLSSEVPQEILPCDPFRPSLLCPQPQALAQSVLQAEIDCDLHSSLQPLL
FPLDTNCSGVFRHRSISPSLLDTAISEEARQGPSLEEEQEQEPLPGSTGRRHTLAEVSTHFSPLNPPC
IIVSSSATASPSEGTSSDCLPFSASEGPAGLGSGLATPGLLGTSSPVRLASPFLGSQSATPVLQTQAGL
GTAVLPPVVSFQEGRRASDTSLTQGLKAFRQQLRKNARTKGFLGLNLIKGLARQVCSSVRTPRGGMSTFH
TPAPSSGLQGCTTSNREGRSLLLEEVHQQRLLQLQHSSSTAAASSGCQGPQLSPVPYVLAPCDLVS
IPLLPTPLLQAGMSPVASAAHLLDTHLHISAGPVALPTGPLPQCLTRLSPGCDPAGLPQGDCEMEDLTSG
QRGTFVLVQ

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/ja1323_g10.zip

Restriction Sites: Sgfl-Mlul

Cloning Scheme:


ACCN: NM_010831

ORF Size: 2337 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_010831.3](#)

RefSeq Size: 4512 bp

RefSeq ORF: 2340 bp

Locus ID: 17691


UniProt ID: [Q60670](#)

Cytogenetics: 17 17.25 cM

MW: 85.6 kDa

Gene Summary: Serine/threonine-protein kinase involved in various processes such as cell cycle regulation, gluconeogenesis and lipogenesis regulation, muscle growth and differentiation and tumor suppression. Phosphorylates HDAC4, HDAC5, PPME1, SREBF1, CRTC1/TORC1 and CRTC2/TORC2. Acts as a tumor suppressor and plays a key role in p53/TP53-dependent anoikis, a type of apoptosis triggered by cell detachment: required for phosphorylation of p53/TP53 in response to loss of adhesion and is able to suppress metastasis. Part of a sodium-sensing signaling network, probably by mediating phosphorylation of PPME1: following increases in intracellular sodium, SIK1 is activated by CaMK1 and phosphorylates PPME1 subunit of protein phosphatase 2A (PP2A), leading to dephosphorylation of sodium/potassium-transporting ATPase ATP1A1 and subsequent increase activity of ATP1A1. Acts as a regulator of muscle cells by phosphorylating and inhibiting class II histone deacetylases HDAC4 and HDAC5, leading to promote expression of MEF2 target genes in myocytes. Also required during cardiomyogenesis by regulating the exit of cardiomyoblasts from the cell cycle via down-regulation of CDKN1C/p57Kip2. Acts as a regulator of hepatic gluconeogenesis by phosphorylating and repressing the CREB-specific coactivators CRTC1/TORC1 and CRTC2/TORC2, leading to inhibit CREB activity. Also regulates hepatic lipogenesis by phosphorylating and inhibiting SREBF1. In concert with CRTC1/TORC1, regulates the light-induced entrainment of the circadian clock by attenuating PER1 induction; represses CREB-mediated transcription of PER1 by phosphorylating and deactivating CRTC1/TORC1.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR219495