

## Product datasheet for **MR219563L4V**

### Mamstr (NM\_172418) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Mamstr (NM_172418) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Mamstr
Synonyms:	2810022D01Rik; 5430432N15Rik; AW743872; Mastr
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_172418
ORF Size:	1263 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR219563).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_172418.2</a> , <a href="#">NP_766006.2</a>
RefSeq Size:	2179 bp
RefSeq ORF:	1266 bp
Locus ID:	74490
UniProt ID:	<a href="#">Q0ZCJ7</a>
Cytogenetics:	7 B3
Gene Summary:	Transcriptional coactivator. Stimulates the transcriptional activity of MEF2C. Stimulates MYOD1 activity in part via MEF2, resulting in an enhancement of skeletal muscle differentiation.[UniProtKB/Swiss-Prot Function]


[View online »](#)