

## Product datasheet for **MR219596**

### Spag11b (NM\_001039563) Mouse Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Spag11b (NM\_001039563) Mouse Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** Spag11b  
**Synonyms:** EG546038; Ep2c/h; Spag11c/h  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**ORF Nucleotide Sequence:** >MR219596 ORF sequence  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGATCGCC**

ATGATACCAAGACTCCTTCCTTTCTTTGTGTCAGTCTACTTTTTGCGGCCCTGCTGTTCCAGGCTTGTCCA  
 ATGCCTCGTCCATTAACCACTTAGTCACTGAACCTCCAGTTTCCCAAGGACGAATTCCTGCTCGAGG  
 AGTAAATGGGTCTCAGTTGTTACATCACAGAGTAAACGATTACCTCCACGCACCCCTTACCACGAG  
 CCTGAACCAAATTACAAATTGCAACTGCAAGAAAAGTGAAGGACAATGTCAAGAATACTGTAATTTCA  
 TGGAAACACAAGTGGGCTACTGTTCAAAAAGAAAGAACCTGCTGCTTACATCCGTTT

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >MR219596 protein sequence  
 Red=Cloning site Green=Tags(s)

MIPRLLPFFVSLLF AALLFPGLSNASSINHLVTEPPSFPKDEF PARGVNGS QLLHHRVKRLPPRTPPYHE  
 PEPNYQIVNCKKSEGQCQEYCNFMETQVGYCSKKKEPCCLHPF

**TR**TRPLEQKLI SEEDLAANDILDYKDDDDKV

**Restriction Sites:** Sgfl-MluI



[View online »](#)

Cloning Scheme:



ACCN: NM\_001039563

ORF Size: 342 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:
1. Centrifuge at 5,000xg for 5min.
  2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
  3. Close the tube and incubate for 10 minutes at room temperature.
  4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
  5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.


RefSeq: [NM\\_001039563.1](#), [NM\\_001039563.2](#), [NM\\_001039563.3](#), [NP\\_001034652.2](#)

RefSeq Size: 642 bp

RefSeq ORF: 342 bp

Locus ID: 546038  
 UniProt ID: [Q8K4N2](#)  
 Cytogenetics: 8 A1.3  
 MW: 12.9 kDa

**Product images:**



Circular map for MR219596