

Product datasheet for **MR219601L4V**

Kif3b (NM_008444) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Kif3b (NM_008444) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Kif3b
Synonyms:	AI854312; AW549267; mKIAA0359
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_008444
ORF Size:	2241 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR219601).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_008444.4 , NP_032470.3
RefSeq Size:	5647 bp
RefSeq ORF:	2244 bp
Locus ID:	16569
UniProt ID:	Q61771
Cytogenetics:	2 75.41 cM
Gene Summary:	Involved in tethering the chromosomes to the spindle pole and in chromosome movement. Microtubule-based anterograde translocator for membranous organelles. Plus end-directed microtubule sliding activity in vitro (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)