

Product datasheet for MR219676

Olfr1535 (NM_207572) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Olfr1535 (NM_207572) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Olfr1535
Synonyms:	MOR256-36; MOR256-63; Olfr1365
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>MR219676 representing NM_207572 Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCGCGATCGCC

ATGAACAGAGCTAATGGCAGCATCTTTTCTGGATTTATTCTTCTGGGCTTCTCCAACAGACCTCAGCTGG
 AAACAGCTCTCTTCTCGCTATGTTGCTCATCTATTTCTTGAGCTTCCTAGGCAACAGCACCATCATATT
 GCTGTGAGTTATAGACCCTCACCTCCACACCCCATGATTTCTCCTTTCTAATCTCTCCTTTATGGAT
 CTTTGTGTTGACAGCTCGCACTGTCCCTCAGACGCTGTTAACTCAAGGGGAAGGACAAGACCATCACTT
 ATGGTGGCTGTGTGACCAACTCTTATTGCCCTGGGACTTGGGGAGTAGAGTGTATACTCTTATCTGT
 GATGGCTTATGATCGCTATGCTGCTGTCTGCCGACCACTCCACTACATGGTAATCATGCACCCACAGCTT
 TGCATAAAGATGGTTATAACTGCTTGTCTACAGGGTTTGCATAAAGATGGTTATAACTGCTTGTCTACAG
 GGTTTGGCAATTCTGTGGTACAGACAGCATTGACCATGAGTCGTCCTCTGTGGCAGAAATCAACTGGA
 CCATTTCTTGTGAAGTTCCAGTGATGCTGAAGCTAACCTGTGCTAACTCCTTCATTAATGAAGCTGAA
 CTCTTTGCTGTCAGTGTCTTCTTGGTGGTACCCCTCTCACTCATATTAGTGTCTATGGCTACATTA
 CTCATGCAGTTCTAAAGATCAAGTCAGCACAGGGAAGGCGGAAGGCCCTTTGGAACCTGTGGCTCTCACCT
 TCTGGTAGTGATTATCTTTTTTGGCACACTCATTTCCATGTACCTCCAGCCTCCCTCCAGTTACTCACAG
 GACATGAACAAAAGCATGGCACTTTTTTATACTCTGGTGACTGCTCTGTTGAATCCCCTGATTTACACAC
 TGAGGAACAGGGAAGTCAAAAGTCATTAAGGAGAATAATGGGTAGAAGGACAGATTGACAAAAGAGT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA


[View online »](#)

Protein Sequence: >MR219676 representing NM_207572
 Red=Cloning site Green=Tags(s)


MNRANGSIFSGFILLGFSNRPQLETALFFAMLLIYFLSFLGNSTIILLVIDPHLHTPMYFFLSNLSFMD
 LCLTARTVPQTLFNFKGDKTITYGGCVTQLFIALGLGGVECILLSVMAYDRYAACRPLHYMIMHPQL
 CIKMVITACSQGLHKDGYNCLLTGFGNSVVQTALMSRPLCGRNQLDHFFCEVPVMLKLTANSFINEAE
 LFAVSVFFLVPLSLILVSYGYITHAVLKIKAQGRRKAFGTGSHLLVVIIFGTLISMYLQPPSSYSQ
 DMNKSMAFYTLVTALLNPLIYTLRNREVKSALRRIMGRRTDSTKS

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mm9072_b12.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_207572

ORF Size: 978 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_207572.1](#), [NP_997455.1](#)

RefSeq Size: 981 bp

RefSeq ORF: 981 bp


Locus ID: 404335

Cytogenetics: 13 A3.1

MW: 36.4 kDa

Gene Summary: Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008]

Product images:



Circular map for MR219676