

Product datasheet for MR219700

Gng12 (NM_001177559) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Gng12 (NM_001177559) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Gng12
Synonyms: 2010305F15Rik; AA536815; AI115529; AI314170; AI842738
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >MR219700 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGTCCAGCAAGACGGCAAGCACCAACAGCATAGCCCAAGCCAGGAGAACTGTGCAGCAGCTGAGATTGG
AAGCCTCCATCGAAAGAATAAAGGTCTCAAAGCATCAGCAGACCTGATGTCATACTGTGAGGAGCATGC
CCGGAGCGACCCCTGCTGATGGGCATACCGACCTCAGAAAACCCGTTCAAGGATAAGAAGACCTGCATC
ATCTTA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR219700 protein sequence
 Red=Cloning site Green=Tags(s)
 MSSKTASTNSIAQARRTVQQLRLEASIERIKVSKASADLMSYCEEHARSDPLLMGIPTSENPFKDKKTCI
 IL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI



[View online »](#)

Cloning Scheme:


ACCN: NM_001177559

ORF Size: 219 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.


RefSeq: [NM_001177559.1](#), [NP_001171030.1](#)

RefSeq Size: 4247 bp

RefSeq ORF: 219 bp

Locus ID: 14701
UniProt ID: [Q9DAS9](#)
Cytogenetics: 6 30.68 cM
MW: 8 kDa
Gene Summary: Guanine nucleotide-binding proteins (G proteins) are involved as a modulator or transducer in various transmembrane signaling systems. The beta and gamma chains are required for the GTPase activity, for replacement of GDP by GTP, and for G protein-effector interaction. [UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR219700