

Product datasheet for MR219730

Fxyd1 (NM_019503) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Fxyd1 (NM_019503) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Fxyd1
Synonyms: 0610012C17Rik; 1110006M24Rik; P; Plm; Pml
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >MR219730 representing NM_019503
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGGGGCAATGGCATCTCCCGCCACATCCTGGCTCTGTGTGTGTCTCTCCATGGCCAGTGCAG
AAGCTCCACAGGAACCGGATCCATTCACCTACGATTACCACACCCTGCGGATCGGCGCCCTCACTATCGC
TGGGATCCTTTCATCTTGGGCATCCTTATCATCCTTAGCAAGAGATGTCGATGCAAATCAACCAACAG
CAGAGAACTGGGAACCCGACGAAGAGGAGGAACTTCCGCAGCTCCATCCGCCGTCTGTATCCCGCA
GGCGG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR219730 representing NM_019503
Red=Cloning site Green=Tags(s)

MGAMASPGHILALCVCLLSMASAEAPQEPDPFTYDYHTLRIGGLTIAGILFILGILIIILSKRCRCKFNQQ
QRTGEPDEEEGTFRSSIRRLSSRRR


TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM_019503

ORF Size: 285 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_019503.4 , NP_062376.2
RefSeq Size:	575 bp
RefSeq ORF:	288 bp
Locus ID:	56188
UniProt ID:	Q9Z239
Cytogenetics:	7 B1
MW:	11 kDa
Gene Summary:	This gene encodes a member of the FXYD family of small membrane proteins that share a 35-amino acid signature sequence domain, beginning with the sequence PFXVD and containing 7 invariant and 6 highly conserved amino acids. The protein encoded by this gene is a plasma membrane substrate for several kinases, including protein kinase A, protein kinase C, NIMA kinase, and myotonic dystrophy kinase. It is thought to form an ion channel or regulate ion channel activity and act as an accessory protein of Na,K-ATPase. Alternatively spliced transcript variants have been described. [provided by RefSeq, Sep 2009]