

Product datasheet for **MR219749L3V**

Hpn (NM_008281) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Hpn (NM_008281) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Hpn
Synonyms:	Hlb32; Hlb320
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_008281
ORF Size:	1248 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR219749).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_008281.4 , NP_032307.2
RefSeq Size:	1770 bp
RefSeq ORF:	1251 bp
Locus ID:	15451
UniProt ID:	O35453
Cytogenetics:	7 B1



[View online »](#)

Gene Summary:

This gene encodes a type II transmembrane serine protease that may function in diverse processes, including regulation of cell growth. Deficiency in this gene results in hearing loss. The protein is cleaved into a catalytic serine protease chain and a non-catalytic scavenger receptor cysteine-rich chain, which associate via a single disulfide bond. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Jan 2013]