

Product datasheet for MR219803

Elavl2 (NM_207685) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Elavl2 (NM_207685) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Elavl2
Synonyms:	Hub; mel-N1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>MR219803 representing NM_207685 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGACAAAGCAATATGAGGTTGCTGTGCAGGTAATTGCTGCCATGGAAACACAACGTCTAATGGGCCAA
CTTGCAATAACACAGCCAATGGTCCAACACTGTAACAACAACGTCTCCTCACCAGTTGACTCCGGGAA
CACAGAGGACAGCAAGACCAACCTAATAGTCAACTACCTCCTCAGAACATGACACAGGAAGAATAAG
AGTCTTTTGGGAGCATTGGCGAGATAGAGTCTGTAAGCTTGTAAAGACAAAAAACAGGGCAGAGCT
TGGGTTACGGATTGTGAACTACATTGACCCCAAGGATGCAGAGAAAGCTATCAACACTCTGAATGGATT
GAGACTTCAAACCAAAACAATAAAAGTTTCTATGCTCGCCCAAGCTCAGCCTCAATCAGAGATGCAAA
TTATACGTGACGGGCTTCCAAAAACATGACCCAGAAGGAGCTGGAACAGCTGTTTTACAGTATGGGC
GCATCATTACTTCTCGGATTCTTGTGACACAGGCTCACTGGCATATCAAGGGGTGTAGGGTTATTCGATT
TGACAAGCGGATTGAGGCAGAAGAAGCTATCAAAGGCCTAAATGGCCAGAAACCTCTGGTGCCACAGAG
CCAATCACTGTAAAGTTTGCTAATAATCCAAGCAAAAAACCAATCAGGCCATCCTTTCCAGCTGTACC
AGTCTCCAACAGAAAGGTACCCAGGACCGCTAGCTCAGCAGGCACAGCGTTTTAGGTTGGACAATCTGCT
CAATATGGCTTATGGAGTGAAGAGGTTTTCTCCAATGACGATCGATGGAATGACCACTCTGGCTGGAATT
AATATTCCTGGCACCCCTGGAACAGGCTGGTGCATATTTGTGTACAACCTGGCGCCTGATGCAGATGAAA
GTATCCTGTGGCAAATGTTTGGGCCTTTTGGAGCTGTACCAACGTGAAGGTATCCGTGACTTTAACAC
CAATAAGTGCAAAGGTTTTGGATTTGTGACTATGACAAACTATGATGAGGCGGCCATGGCGATAGCTAGC
CTCAATGGATACCGCTGGGAGACAGAGTCTGCAAGTCTCCTTTAAGACAAAACAAACGCACAAAGCC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >MR219803 representing NM_207685
 Red=Cloning site Green=Tags(s)

MTKQYEVAVQVIAAMETQLSNGPTCNNTANGPTTVNNNCSSPVDSGNTEDSKTNLIVNYLPQNMTQEELK
 SLFGSIGEIESCKLVRDKITGQSLGYGFVNYIDPKDAEKAINLNLRLQTKTIKVSYPSSASIRDAN
 LYVSGLPKMTQKELEQLFSQYGRITTSRILVDQVTGISRGVGFIKRFDRKIEAEEAIKGLNGQKPPGATE
 PITVKFANNPSSQKTNQAILSQLYQSPNRRYPGPLAQQAQRFRLDNLLNMAYGVKRFSPMTIDGMTSLAGI
 NIPGHPGTGWCIFVYNLAPDADESILWQMFPGFVAVTNVAVKIRDFNTNCKKGFVVTMTNYDEAAMAIAS
 LNGYRLGDRVLQVSFKTNKTHKA

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF

ACCN: NM_207685

ORF Size: 1119 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_207685.2](#)

RefSeq Size: 3889 bp

RefSeq ORF: 1122 bp


Locus ID: 15569

Cytogenetics: 4 42.36 cM

MW: 41.5 kDa

Gene Summary: Binds RNA. Seems to recognize a GAAA motif. Can bind to its own 3' UTR, the FOS 3' UTR and the ID 3' UTR.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR219803