

Product datasheet for **MR219825L3V**

Rpf2 (NM_001042556) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Rpf2 (NM_001042556) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Rpf2
Synonyms:	2810470K21Rik; AU040229; AU043242; Bxdc1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001042556
ORF Size:	819 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR219825).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001042556.1 , NP_001036021.1
RefSeq Size:	1830 bp
RefSeq ORF:	822 bp
Locus ID:	67239
UniProt ID:	Q9J180
Cytogenetics:	10 B1
Gene Summary:	Involved in ribosomal large subunit assembly. May regulate the localization of the 5S RNP/5S ribonucleoprotein particle to the nucleolus.[UniProtKB/Swiss-Prot Function]



[View online »](#)