

Product datasheet for **MR219831L3V**

Fancb (NM_001146081) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Fancb (NM_001146081) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Fancb
Synonyms:	BC022692
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001146081
ORF Size:	2559 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR219831).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001146081.1 , NP_001139553.1
RefSeq Size:	2984 bp
RefSeq ORF:	2562 bp
Locus ID:	237211
UniProt ID:	Q5XJY6
Cytogenetics:	X F5
Gene Summary:	DNA repair protein required for FANCD2 ubiquitination.[UniProtKB/Swiss-Prot Function]



[View online »](#)