

Product datasheet for MR219834L3

Ano10 (NM_133979) Mouse Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Ano10 (NM_133979) Mouse Tagged Lenti ORF Clone
Tag:	Myc-DDK
Symbol:	Ano10
Synonyms:	AI604832; Tmem16k
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR219834).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF.


ACCN:	NM_133979
ORF Size:	1977 bp



[View online »](#)

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_133979.2 , NP_598740.1
RefSeq Size:	2653 bp
RefSeq ORF:	1980 bp
Locus ID:	102566
UniProt ID:	Q8BH79
Cytogenetics:	9 F4
Gene Summary:	This gene encodes a member of the anoctamin family, which in mammals is comprised of 10 members. Anoctamin proteins are proposed to have eight transmembrane domains with both termini facing the cytoplasm and a C-terminal domain of unknown function. While some members have been characterized as calcium-activated chloride channels, this protein is reported to inhibit anion conductance. Alternative splicing results in multiple transcript variants that encode different protein isoforms. [provided by RefSeq, Dec 2012]

Product images:



Circular map for MR219834L3