

Product datasheet for **MR219853L4V**

IL3ra (NM_008369) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	IL3ra (NM_008369) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	IL3ra
Synonyms:	CD123; CDw123; SUT-1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_008369
ORF Size:	1191 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR219853).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_008369.1 , NP_032395.1
RefSeq Size:	1394 bp
RefSeq ORF:	1191 bp
Locus ID:	16188
UniProt ID:	P26952
Cytogenetics:	14 7.08 cM
Gene Summary:	In mouse, there are two classes of high-affinity IL3 receptors. One contains an IL3-specific beta subunit and the other contains the beta subunit also shared by high-affinity IL5 and GM-CSF receptors.[UniProtKB/Swiss-Prot Function]



[View online »](#)