

Product datasheet for **MR219870L4V**

Cep68 (NM_172260) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Cep68 (NM_172260) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Cep68
Synonyms:	6030463E10Rik; AI481761; BC027174; Kiaa0582
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_172260
ORF Size:	2199 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR219870).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_172260.3 , NP_758464.2
RefSeq Size:	4810 bp
RefSeq ORF:	2202 bp
Locus ID:	216543
UniProt ID:	Q8C0D9
Cytogenetics:	11 12.92 cM
Gene Summary:	Involved in maintenance of centrosome cohesion, probably as part of a linker structure which prevents centrosome splitting. Required for localization of CDK5RAP2 to the centrosome during interphase.[UniProtKB/Swiss-Prot Function]



[View online »](#)