

Product datasheet for **MR219872**

Krt36 (NM_001174099) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Krt36 (NM_001174099) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Krt36
Synonyms:	HRa-1; K36; Krt1-5; Krt1-22; MHRa-1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>MR219872 representing NM_001174099
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**GCGATCGCC**

ATGGCCACCCAGATCTGCACCCGACCTTCTCTGCTGGGTCTGCCAAGGCCTCTGTGGCACTTCAGGGA
 GCTTCTCTCGGATATCTCCATCCATTCCATGGGCTCCTGCAGGGCCCCAGTTGGTTGGCACTGTGGG
 GTCGCTGTCTCTCCTCAGGACAGGCTTCTCTGGCCCTGGGAGCTGCTTGCAGGGCTCCTACCTGTCTCT
 GGGTGCCACTCTCTGGCTTTGCTGGGAGCGGGGATGGTTCTGCAGGGCGCCTTCAATGGCAACGAGA
 AGGCGACCATGCAGATCTTGAATGACCGACTGGCCAACTACCTGGAGAAGGTGCGGCAGCTGGAGCAGGA
 GAACACACAGCTGGAGTGCCGCATCCGAGAGTGGTACGAGTGTGAGATCCCATACATCTGCCCGACTAC
 CAGTCTACTTCAAGACTGCTGAGGAGTGCAGCAGAAGATCCTGCTGACCAAGTCTGAGAATGCCAGGC
 TGATATTACAAATTGACAATGCCAAGCTGGCTGCTGATGACTCCGGACCAAGTATGAGACGGAGTTGTC
 CCTGCGACAGCTGGTAGAGGCTGACATCAATGGTCTACGCCGGATCTTGGATGAGCTGACCTTGTAAG
 GCTGACCTGGAGGCACAGGTTGAGTCCCTGAAGGAGGAGCTGCTGTCTCAAGAGGAACCATGAAGAGG
 AAGTCAATGCCCTCCGTTCCAGCTTGGGACCGACTGAATGTGGAGGTGGATGCTGCTCCCCAGTGGA
 TCTCAACAAGATACTGGATGACATGCGATGCCAGTATGAGACGCTGGTGGAGAATAACCGTAGAGATGTG
 GAGGCTTGGTTCAACACTCAGACCGAGGAGCTGAACCAGCAGGTGGTATCCAGCTCAGAGCAGCTGCAGT
 GTTGCCAGACAGAGATCATTGAGCTGAGACGCACAGTCAACTCCCTGGAGATCGAGCTGCAGGCTCAACA
 GAGCATGAGGAATTCTCTGGAATCTACCCTGGCTGAGACAGAGGCCCGCTACAGCTCCAGCTGGGCCAG
 ATGAGTATCAGGTGTTACTGGATGTCAAGGCCCGGCTGGAGTCAAGATGCCCACCTACCGTCCGCTGCTGCT
 GGATGGGAGGACTGCAAGCTTCTGCCACCCTTGTCCACGGAATGCAAGCCTGCTGTAAGAGTACCT
 TATATCCCTTCTACAACCTGCACTCCAGCTGGGCCCTGCACCCAGCTGGGCCCTGCACCCAGCTCCTC
 AAGTCAGCACCCAGATCCGGACCATCACCGAGGAGATCAGGGATGGAAGAGTCATCTCTCCAGGGAGCA
 CGTGGTGCCCGGCAATG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>MR219872 representing NM_001174099
 Red=Cloning site Green=Tags(s)

MATQICTPTFSAGSAKGLCGTSGSFSRISSIHSMGSCRAPSLVGTVGSVSSFRITGFSGPGSCLPGSYLSS
 GCHSSGFAGSGGWFCEGAFNGNEKATMQILNDRLANYLEKVRQLEQENTQLECRIREWYECQIPYICPDY
 QSYFKTAEELQKILLTKSENARLILQIDNAKLAADDFRTKYETELSLRQLVEADINGLRRILDELTLCK
 ADLEAQVESLKEELLCLKRNHEEEVNALRSQGDRLNVEVDAAPPVDLNKILDDMRCQYETLVENNRDV
 EAWFNTQTEELNQQVSSSEQLQCCQTEIIELRRTVNSLEIELQAQQSMRNSLESTLAETEARYSSQLGQ
 MQCLITNVESQLAEIRCDLERQNHEYQVLLDVKARLESEIATYRRLLDGEDCKLPAHPCSTECKPAVRVP
 YIPSTTCTPAGPCTPAGPCTPAPQVSTQIRTITEEIRDGRVISSREHVVPAM

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms:

https://cdn.origene.com/chromatograms/mm9075_f02.zip

Restriction Sites:

Sgfl-Mlul

Cloning Scheme:


ACCN: NM_001174099

ORF Size: 1419 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001174099.1](#), [NP_001167570.1](#)

RefSeq Size: 1666 bp

RefSeq ORF: 1422 bp


Locus ID: 16673

UniProt ID: [B1AQ75](#)

Cytogenetics: 11 D

MW: 52.8 kDa

Product images:



Circular map for MR219872