

Product datasheet for **MR219883A1V**

Mouse Ccny (NM_026484) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Mouse Ccny (NM_026484) AAV Particle
Tag:	Myc-DDK
Symbol:	Ccny
Synonyms:	1700025H17Rik; 3110050L10Rik; 4631402G10Rik; 5730405I09Rik
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>MR219883 representing NM_026484 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCCGCC**CGGATCGCC**

ATGGGGAACACAACCTTCGTGCTGCGTGTCCAGCCCCAAGCTCCGGAGGAATGCGCACTCGGGCTGGAGTCCCTACCGCCGGACCGGACCTCAGCCGTGAGGACTGGCTGCAATCTGCAGCACATCAGCGACCGGGAAAACATCGACGATTTAAACATGGAATTCATCCTTCAGACCATCCTCGGGCCAGCACAATATTCCTCAGCAAATCTCAGACAGACGTGAGAGAAAAACGCAAGAGTCTCTTCATTAACCATCATCCTCCAGGGCAGACATCAAGGAAATACAGCTCTTGCTCCACCATCTTCCTAGATGACAGCACAGTCAGCCAGCCAAACCTCAAATGTACACCATCAAATGTGTAGCTCTTGCGATATATTACCACATCAAAAACAGGGACCCAGATGGAAGGATGCTCTTAGATATTTTTGATGAAAACCTCCACCCTCTTTCGAAATCCGAAGTGCCACCTGATTATGACAAACACAACCCAGAACAGAAGCAAATTTACCGGTTTCGTTTCGGACTCTTCAGTGTGCTCAGCTGACAGCTGATGTGCCATTGTCACCCTGGTATACCTTGAAAGGCTCTTGACATACGCAGAGATTGATATCTGTCCAGCCAATTGGAAGCGCATTGTTCTAGGGCAATCCTGCTGGCCTCTAAGGTGTGGATGACCAAGCTGTGTGGAATGTGGATTACTGCCAGATCCTGAAAGACATCACAGTGGAGGACATGAATGAACTGGAGAGGCAGTTTCTGGAATTGCTACAGTTAACATCAATGTTCCCTTCAGTGTGTATGCCAAGTATTATTTTGTCTTCGTTCTCTGGCAGAAGCAAACAACCTGAGCTTTCTCTGGAGCCCCTAAGCAGAGAGCGTGCCACAAGCTTGAGGCCATCTCTCGTCTCTGTGAGGACAAGTACAAGGACCTGAGGAAACCCATGAGGAAGCGTTCAGCAAGTGTGACAACCTGATTCTGCCAAGGTGGTCTCCAGCCATCATCTCC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATTACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>MR219883 representing NM_026484 Red=Cloning site Green=Tags(s) MGNTTSCCVSSPKLRRNAHSRLESYRPDTDL SREDTGCNLQHISDRENIDDLNMEFNPSDHPRASTIFL SKSQTDVREKRKSLFINHHPPGQTSRKYSSCSTIFLDDSTVSQP NLKYTIKCV ALAIYYHIK NRDPDGRM LLDIFDENLHPLSKSEVPPDYDKHNPEQKQIYRFVRTLF SAAQLTAECAIVTLVYLERLLTYAEIDICPA NWKRIVLGAILLASKVWDDQAVWNVDYCQILKDITVEDMNELE RQFLELLQFNINVPSSVYAKYYFDLRS LAEANNLSFPLEPLSRERAHKLEAISRLCEDKYKDLRKP MRKRSASADNLILPRWSPA IIS TRTRPLEQKLI SEEDLAANDILDYKDDDDKV
Species:	Mouse
Serotype:	AAV-2
ACCN:	NM_026484
ORF Size:	1026 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_026484.3
RefSeq Size:	4017 bp
RefSeq ORF:	1026 bp
Locus ID:	67974
UniProt ID:	Q8BGU5
Cytogenetics:	18 A1
MW:	39.4 kDa