

## Product datasheet for **MR219883L1V**

### Ccny (NM\_026484) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Ccny (NM_026484) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Ccny
Synonyms:	1700025H17Rik; 3110050L10Rik; 4631402G10Rik; 5730405I09Rik
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_026484
ORF Size:	1026 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR219883).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_026484.3</a>
RefSeq Size:	4017 bp
RefSeq ORF:	1026 bp
Locus ID:	67974
UniProt ID:	<a href="#">Q8BGU5</a>
Cytogenetics:	18 A1



[View online »](#)

**Gene Summary:**

Positive regulatory subunit of the cyclin-dependent kinase CDK14/PFTK1. Acts as a cell-cycle regulator of Wnt signaling pathway during G2/M phase by recruiting CDK14/PFTK1 to the plasma membrane and promoting phosphorylation of LRP6, leading to the activation of the Wnt signaling pathway (By similarity). Recruits CDK16 to the plasma membrane (By similarity). Positive regulatory subunit of the cyclin-dependent kinase CDK16.[UniProtKB/Swiss-Prot Function]