

Product datasheet for **MR219888L3V**

Anks1b (NM_001177398) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Anks1b (NM_001177398) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Anks1b
Synonyms:	AIDA-1b; C030032C09Rik; E530015N03Rik; ENSMUSG00000078431; Gm1555; Gm10937
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001177398
ORF Size:	1353 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR219888).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001177398.1
RefSeq Size:	4728 bp
RefSeq ORF:	1356 bp
Locus ID:	77531
UniProt ID:	Q8BIZ1
Cytogenetics:	10 C2
Gene Summary:	Isoform 2 may participate in the regulation of nucleoplasmic coilin protein interactions in neuronal and transformed cells.[UniProtKB/Swiss-Prot Function]



[View online »](#)