

## Product datasheet for **MR219891L3V**

### Akirin2 (NM\_001007589) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Akirin2 (NM_001007589) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Akirin2
Synonyms:	2700059D21Rik; AA114675; AA522011; Akirin-2; AU019887
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001007589
ORF Size:	603 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR219891).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001007589.3</a> , <a href="#">NP_001007590.2</a>
RefSeq Size:	1729 bp
RefSeq ORF:	606 bp
Locus ID:	433693
UniProt ID:	<a href="#">B1AXD8</a>
Cytogenetics:	4 A5



[View online »](#)

**Gene Summary:**

Forms a complex with YWHAB that acts to repress transcription of DUSP1 (By similarity). Required for embryonic development and the innate immune response. Downstream effector of the Toll-like receptor (TLR), TNF and IL-1 beta signaling pathways leading to the production of IL-6.[UniProtKB/Swiss-Prot Function]