

## Product datasheet for **MR219914L3V**

### **Cacng5 (NM\_080644) Mouse Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Cacng5 (NM_080644) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Cacng5
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_080644
ORF Size:	828 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR219914).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_080644.3</a> , <a href="#">NP_542375.1</a>
RefSeq Size:	3678 bp
RefSeq ORF:	828 bp
Locus ID:	140723
UniProt ID:	<a href="#">Q8VHW4</a>
Cytogenetics:	11 E1



[View online »](#)

**Gene Summary:**

Regulates the gating properties of AMPA-selective glutamate receptors (AMPA receptors). Modulates their gating properties by accelerating their rates of activation, deactivation and desensitization. Displays subunit-specific AMPA receptor regulation. Shows specificity for GRIA1, GRIA4 and the long isoform of GRIA2. Thought to stabilize the calcium channel in an inactivated (closed) state (By similarity).[UniProtKB/Swiss-Prot Function]