

Product datasheet for **MR219936L4V**

Parvg (NM_001162500) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Parvg (NM_001162500) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Parvg
Synonyms:	AI413459
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001162500
ORF Size:	1152 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR219936).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001162500.1 , NP_001155972.1
RefSeq Size:	2700 bp
RefSeq ORF:	1155 bp
Locus ID:	64099
UniProt ID:	Q9ERD8
Cytogenetics:	15 E2
Gene Summary:	Probably plays a role in the regulation of cell adhesion and cytoskeleton organization. [UniProtKB/Swiss-Prot Function]



[View online »](#)