

## Product datasheet for **MR219938L4V**


### **Pgs1 (NM\_133757) Mouse Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Pgs1 (NM_133757) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Pgs1
Synonyms:	2610019F11Rik; 4933424M23Rik; SAF
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_133757
ORF Size:	1659 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR219938).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_133757.2</a> , <a href="#">NP_598518.1</a>
RefSeq Size:	2298 bp
RefSeq ORF:	1662 bp
Locus ID:	74451
UniProt ID:	<a href="#">Q8BHF7</a>
Cytogenetics:	11 E2
Gene Summary:	Functions in the biosynthesis of the anionic phospholipids phosphatidylglycerol and cardiolipin.[UniProtKB/Swiss-Prot Function]



[View online »](#)

**Product images:**

[MR219938L4] was used to prepare Lentiviral particles using [TR30037] packaging kit. HEK293T cells were transduced with MR219938L4V particle to overexpress human Pgs1-mGFP fusion protein.