

## Product datasheet for MR219953

### Tssk6 (NM\_032004) Mouse Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Tssk6 (NM_032004) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Tssk6
Synonyms:	Sstk; TSK-6; TSSK-6
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>MR219953 representing NM_032004 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGTCGGGCGACAACTCCTGAGCGAACTTGGCTATAAGCTGGGCCGACGATAGGTGAGGGCAGCTACT  
CCAAGGTGAAAGTGCCACCTCCAAGAAGTATAAAGGGACGGTGGCCATCAAGGTGGTGGACCGCCGGCG  
AGCGCCTCCGGACTTCGTCAACAAGTTCCTGCCTCGAGAGCTGTCCATCCTCGGGGAGTACGGCACCCG  
CACATCGTACACGTTTTTCGAGTTCATCGAAGTGTCAATGGGAAGCTGTACATCGTGATGGAGGCAGCAG  
CCACCGACTGTGCAGGCAGTGCAGCGCAACGGGCGCATCCCGGGTTCGAGGCGCGCAACTCTTCTC  
GCAGATCGCGGGCGCCGTGCGCTACCTGCACGACCACCCTGGTGCACCGCGACCTCAAGTGCGAAAAC  
GTGCTGCTAAGCCCTGACGAGCGCCGCTCAAGATCACGGATTTTCGGTTTCGGCCGCCAGGCGCACGGCT  
ACCCAGACCTGAGCACCCTACTGCGGCTCGGCCGCTATGCGTACCTGAGGTGCTCCTGGGCATCCC  
CTACGACCCTAAGAAATACGACGTGTGGAGCCTGGGCGTTGTGCTCTACGTATGGTACCCGGGTGCATG  
CCCTTCGACGACTCGGACATTGCTGGCTACCCCGGCCAGAAAGCGCGGCGTCTACCTGACGGCC  
TCGAACTGTCGGAACGATGCAAATCCTTGATCGCCGAGCTGTACAGTTCAGCCCGTCTGCCCGTCCCTC  
GGCGGGCCAGGTGGCGCGCAATGGCTGGCTGCGGGCCGGGACTCCGGC

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

**Protein Sequence:** >MR219953 representing NM\_032004  
 Red=Cloning site Green=Tags(s)

MSGDKLLSELGYKLGRTIGEGSYSKVKVATSKKYKGTVAIKVVDRRRAPPDFVNKFLPRELSILRGVRHP  
 HIVHVFEFIEVCNGKLYIVMEAAATDLLQAVQRNGRIPGSQARELFSQIAGAVRYLHDHHLVHRDLKCEN  
 VLLSPDERRVKITDFGFRQAHGYPDLSSTYCGSAAAYASPEVLLGIPYDPKKYDWSLGVVLYVMVTGCM  
 PFDDSDIAGLPRRQKRGVLYPDGLELSERCKSLIAELLQFSPSARPSAGQVARNGWLRAGDSG

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

**Chromatograms:** [https://cdn.origene.com/chromatograms/mm9010\\_e02.zip](https://cdn.origene.com/chromatograms/mm9010_e02.zip)

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**

Cloning sites used for ORF Shuttling:



\* The last codon before the Stop codon of the ORF

**ACCN:** NM\_032004

**ORF Size:** 819 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_032004.1](#), [NP\\_114393.1](#)

**RefSeq Size:** 1327 bp

**RefSeq ORF:** 822 bp

**Locus ID:** 83984


**UniProt ID:** [Q925K9](#)

**Cytogenetics:** 8 B3.3

**MW:** 30.7 kDa

**Gene Summary:** Required for sperm production and function. Plays a role in DNA condensation during postmeiotic chromatin remodeling.[UniProtKB/Swiss-Prot Function]

## Product images:



Circular map for MR219953