

Product datasheet for **MR219970L3V**

Angptl2 (NM_011923) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Angptl2 (NM_011923) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Angptl2
Synonyms:	AI593246; Arp2; AW260363
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_011923
ORF Size:	1479 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR219970).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_011923.4 , NP_036053.2
RefSeq Size:	3444 bp
RefSeq ORF:	1482 bp
Locus ID:	26360
UniProt ID:	Q9R045
Cytogenetics:	2 B
Gene Summary:	Induces sprouting in endothelial cells through an autocrine and paracrine action. [UniProtKB/Swiss-Prot Function]



[View online »](#)