

Product datasheet for **MR219984L3V**

Cdk15 (NM_001033373) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Cdk15 (NM_001033373) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Cdk15
Synonyms:	Als2cr7; Pftk2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001033373
ORF Size:	1302 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR219984).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001033373.2
RefSeq Size:	1477 bp
RefSeq ORF:	1302 bp
Locus ID:	271697
UniProt ID:	Q3V3A1
Cytogenetics:	1 C1.3- C2
Gene Summary:	Serine/threonine-protein kinase that acts like an antiapoptotic protein that counters TRAIL/TNFSF10-induced apoptosis by inducing phosphorylation of BIRC5 at 'Thr-34'. [UniProtKB/Swiss-Prot Function]



[View online »](#)