

Product datasheet for **MR220102L3V**

Serpinf2 (NM_008878) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Serpinf2 (NM_008878) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Serpinf2
Synonyms:	A1747498; Pli; Serpimf2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_008878
ORF Size:	1473 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR220102).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_008878.2 , NP_032904.1
RefSeq Size:	2199 bp
RefSeq ORF:	1476 bp
Locus ID:	18816
UniProt ID:	Q61247
Cytogenetics:	11 B5
Gene Summary:	Serine protease inhibitor. The major targets of this inhibitor are plasmin and trypsin, but it also inactivates matriptase-3/TMPRSS7 and chymotrypsin (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)