

Product datasheet for **MR220156L3V**

Arhgef10l (NM_001112722) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Arhgef10l (NM_001112722) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Arhgef10l
Synonyms:	2810441C07Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001112722
ORF Size:	3723 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR220156).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001112722.1 , NP_001106193.1
RefSeq Size:	4450 bp
RefSeq ORF:	3726 bp
Locus ID:	72754
UniProt ID:	A2AWP8
Cytogenetics:	4 D3
Gene Summary:	Acts as guanine nucleotide exchange factor (GEF) for RHOA, RHOB and RHOC. [UniProtKB/Swiss-Prot Function]



[View online »](#)