

## Product datasheet for **MR220159L4V**

### Eps8l1 (NM\_026146) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Eps8l1 (NM_026146) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Eps8l1
Synonyms:	2310051G19Rik; 4632407K17Rik; AW060268; DRC3; EPS8R1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_026146
ORF Size:	2148 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR220159).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_026146.4</a> , <a href="#">NP_080422.1</a>
RefSeq Size:	3774 bp
RefSeq ORF:	2151 bp
Locus ID:	67425
UniProt ID:	<a href="#">Q8R5F8</a>
Cytogenetics:	7 A1
Gene Summary:	Stimulates guanine exchange activity of SOS1. May play a role in membrane ruffling and remodeling of the actin cytoskeleton (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)