

## Product datasheet for **MR220189L3V**

### Slc35f4 (NM\_029238) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Slc35f4 (NM_029238) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Slc35f4
Synonyms:	4930550L21Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_029238
ORF Size:	1455 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR220189).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_029238.2</a> , <a href="#">NP_083514.1</a>
RefSeq Size:	1881 bp
RefSeq ORF:	1458 bp
Locus ID:	75288
UniProt ID:	<a href="#">Q8BZK4</a>
Cytogenetics:	14 C1
Gene Summary:	Putative solute transporter.[UniProtKB/Swiss-Prot Function]



[View online »](#)