

## Product datasheet for **MR220203**

### **Pja2 (NM\_001025309) Mouse Tagged ORF Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	Pja2 (NM_001025309) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Pja2
Synonyms:	AI447901; AL022700; mKIAA0438; Neurodap1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**ORF Nucleotide Sequence:**

>MR220203 representing NM\_001025309  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGGATCGCC**

ATGTCACAGTACACTGAAAAAGAGCCATCAGTAATGGATCAAGAATCCAGCAAGGCTGCTTGGCCTAAAC  
 CAGCAGGAGGATATCAGACTATCACAGGCAGGAGATATGGAAGAAGACATGCGTATGTCAGTTTTAAGCC  
 ATGTATGACCAGGCATGAAAGAAGCTTGGGTCGGGCTGGTGTACTATGAAGTGTGGAAGTACAGCAGT  
 GTTCCAAAGGAGAACACCTCAGGTTCCAGTTCATTGGATCAAGTCCATCCTGCTTTACCAAATGAACCTA  
 CAGTTGAAAAAGTGAACAGAAATTTCTACTTGGTCCAGCACTGAATCAAAGCACTGAGAGCAGTCC  
 ATCCATTGCCACAGTATGTCATAGTGAGGAAGTCAGGGAGACCTAGAAAAGCAATACAAACCTTCATAAT  
 CGCACTGAGACAGAGCATACTCCAGCAGTTTGAATGTCTCAAGTGTCCAGAATGGAATTATGTTGGTTC  
 ATACTGACTCTTATGATCCAGACAGCAAGCATGATGAGAATGGCTCTTCAACTGGCGCTGAAGCTGT  
 GGAAGGTGGTAGACATCAGAAGGGGTTAGGCCGTGCAGTCTTTGAGCTGGAAAATGGAGAGGCAGAGATC  
 TACGCTGACCTGTCAACATCAGTCCCCTCTCTCAACGGTAAAATAGTGAGGCGTTTTGAAGAGCTAGATT  
 CAGCACCTTTAGAGAAAAGTTCTACTGCCGATGCTGAGCTTGCCATCAGAATGGGCAGGAATTTACAGAG  
 GTCTTCTGAAGATGGCGTTGTTAGAAAAGAGGCGACAAGACGATACTGATCAGGGAAGACAGACAGAAAAT  
 TCAACTGAAGATGCAGACTGTGCTCCAGGCCATGTTGAGCAAAAATACTAGTGATAGAGCCAAACCACATG  
 GAAGTTCTCCTGAACAGGTAGTAAGGCCAAAAGTTAGAAAAGTGATAAGTTCAAGCCAGGTGGACCAAGA  
 AATTGGTTTTAATAGGCACGAGGCTAAGCAGAGAAGCGTTCAGAGATGGAGAGAGGCTTTGGAAGTGGAG  
 GAGTGCAGTTCAGATGACCCTATAATAAGTGTGACGATTATGATGGAGACCACGACTGCATGTTCTTAA  
 CTCCAACCTTACTCAAGAGTTACACAAAGGGAGACAGAACGGAATCGTGTAAATCAGAAAATGGAGCTAC  
 AGCTTCAGGACGGCAAGAATCTCGGGACAATGCCTTTTTGGAATGCCTGCGGAGAGTATTACCAACTCTTT  
 GACAAAGATGAAGACAGTTCAGAGTGCAGTGTGGGAATGGTCTGCTTCTGCTCACCATTCTGCTG  
 GCACAGAAAAAGACCAGTCTCAAGCGATGACAGCTGGGAACTCTGCCAGGAAAGGATGAGAATGACCC  
 TGAGCTGCAGAGTGATAGTAGTGGCCCCGAGGAAGAAAACCAAGAGTTGTCTTTCAGGAGGGGGAACAG  
 ACATCCTTGAAGAGGGAGAGATCCCTTGGTTACAGTACAATGAGGTCAATGAGAGCAGCAGCGATGAAG  
 GGAATGAGCCTGCCAATGAGTTGCACAGCCGGAAGCTTTCATGTTGGATGGGAACAACAACCTTGAGGA  
 CGACTCGAGCGTGAAGACTTGGATGTGGACTGGAGCCTGTTTATGTTTGCAGACGGACTTGA  
 GTTGCTGAAGCTATTTCTTATGTGGATCCTCAGTTCCTTACCTACATGGCACTAGAAGAGCGCTTAGCCC  
 AGGCTATGGAGACTGCGCTGGCACATTTAGAGTCCCTTGCAGTAGACGTTGAGGTGGCCAATCCACCTGC  
 CAGTAAGGAAAGCATTGATGGGCTTCCAGAGACTCTGGTTCTAGAGGACCACACCGCCATCGGCCAGGAG  
 CAGTGTGCTCCCATCTGCTGCAGCGAGTACATCAAGGATGACATTGCCACAGAGCTGCCTTGCCATCACT  
 TCTTTCAAGCCCTGTGTCTCCATCTGGCTACAGAAGTCCGGAACATGCCCTGTGTGCCCGCTCACTT  
 CCCACCTGCAGTATTGACGCATCTGCAGCTGCTTCTGACCTGACCCTGATGCCTCGCTGCAAAAT  
 GACAATGCTGAAGAAGCCCCA

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >MR220203 representing NM\_001025309  
Red=Cloning site Green=Tags(s)


MSQYTEKEPSVMDQESSKAAWPKPAGGYQTITGRRYGRRHAYVSFKPCMTRHERSLGRAGDDYEVLLEDD  
 VPKENTSGSSSLDQVHPALPNEPTVEKSETEISTCGPALNQSTESSPSIATVCHSEEVRETLESNTNLHN  
 RTETEHTPAVCNVSSVQNGIMLVHTDSYDPDSKHDENGLQLGAEAVEGGRHQKGLGRAVFELENGEAEI  
 YADLSPSVPSLNGEISEAFEELDSAPLEKSSTADAELVHQNGQEFQRSSSEDGVVRKRRQDDTDQGRQTEN  
 STEDADCAPGHVEQNTSDRANHGGSSPEQVVRPKVRKVISSSQVDQIEGFNRHEAKQRSVQRWREALVE  
 ECSSDDPIIKCDDYDGDHDCMFLPTYSRVVTQRETERNRVTSENGATASGRQESRDNAFWNACGEYYQLF  
 DKEDSSECSGDGEWSASLPHRFSGTEKDQSSDDSWETLPGKDENDPELQSDSSGPPEENQELSLQEGEQ  
 TSLEEIGEIPWLQYNEVNESSDEGNEPANEFAPQPEAFMLDGNLLEDDSSVSEDLVDVSLFDGADGLG  
 VAEAISYVDPQFLTYMALEERLAQAMETALAHLESLAVDVEVANPPASKESIDGLPETLVLEDHTAIGQE  
 QCCPICSEYIKDDIATELPCHHFFHKPCVSIWLQKSGTGPVCRRHFPFAVIDASAAASSDPDPDASPAN  
 DNAAEAP

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

**Restriction Sites:**

SgfI-MluI

**Cloning Scheme:**



**ACCN:** NM\_001025309

**ORF Size:** 2121 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001025309.1](#), [NP\\_001020480.1](#)

**RefSeq Size:** 4706 bp

**RefSeq ORF:** 2124 bp

**Locus ID:** 224938


**UniProt ID:** [Q80U04](#)

**Cytogenetics:** 17 E1.1

**MW:** 78.4 kDa

**Gene Summary:** Has E2-dependent E3 ubiquitin-protein ligase activity. Responsible for ubiquitination of cAMP-dependent protein kinase type I and type II-alpha/beta regulatory subunits and for targeting them for proteasomal degradation. Essential for PKA-mediated long-term memory processes. Through the ubiquitination of MFHAS1, positively regulates the TLR2 signaling pathway that leads to the activation of the downstream p38 and JNK MAP kinases and promotes the polarization of macrophages toward the pro-inflammatory M1 phenotype.[UniProtKB/Swiss-Prot Function]

**Product images:**



Circular map for MR220203