

Product datasheet for **MR220211L3V**

Mthfs (NM_026829) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Mthfs (NM_026829) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Mthfs
Synonyms:	1110034I12Rik; 2310020H23Rik; AI119695
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_026829
ORF Size:	609 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR220211).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_026829.2 , NP_081105.1
RefSeq Size:	768 bp
RefSeq ORF:	612 bp
Locus ID:	107885
UniProt ID:	Q9D110
Cytogenetics:	9 E3.1



[View online »](#)

Gene Summary:

Contributes to tetrahydrofolate metabolism. Helps regulate carbon flow through the folate-dependent one-carbon metabolic network that supplies carbon for the biosynthesis of purines, thymidine and amino acids. Catalyzes the irreversible conversion of 5-formyltetrahydrofolate (5-CHO-H(4)PteGlu) to yield 5,10-methenyltetrahydrofolate (By similarity).[UniProtKB/Swiss-Prot Function]