

Product datasheet for **MR220245L4V**

Pla2g4f (NM_001024145) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Pla2g4f (NM_001024145) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Pla2g4f
Synonyms:	4732472I07Rik; Gm647; Pla2zeta
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001024145
ORF Size:	2565 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR220245).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001024145.2
RefSeq Size:	3125 bp
RefSeq ORF:	2568 bp
Locus ID:	271844
UniProt ID:	Q50L41
Cytogenetics:	2 E5
Gene Summary:	Calcium-dependent phospholipase A2 that selectively hydrolyzes glycerophospholipids in the sn-2 position. Has higher enzyme activity for phosphatidylethanolamine than phosphatidylcholine.[UniProtKB/Swiss-Prot Function]



[View online »](#)