

Product datasheet for MR220254

S100a7a (NM_199422) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: S100a7a (NM_199422) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: S100a7a
Synonyms: AY465109; Gm1020; S100A7f; S100A7L1; S100a15; S100a15a; S100a1711
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >MR220254 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCCGCATCGCC

ATGCCAGACACACCAGTGGAGGACTCCCTCTTCCAAATCATACTGCTTCCATCACTATGCTGCCCGGG
 AAGGGACAAGGAGACCTTGTCCCTGGAGGAGTTGAAAGCTCTGCTTTGGATAGTGTGCCTCGCTTCAT
 GGACACCTTGGCCCGCAGGCAGCCATACTACATCACAGAGCTGTCCGGGCAGCTGACAAAACAAGGAC
 AACCAGATCTGCTTTGATGAGTTCCTGTACATCTGGGCAAGCTGGTGAAGGACTACCATCTCCAGTTCC
 ACCGGCAGTTGTGTGCACACTACTGTACTGAGCACAGCCTCTAC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR220254 protein sequence
 Red=Cloning site Green=Tags(s)
 MPDTPVEDSLFQIIHCFHHYAAREGDKETLSLEELKALLLDSVPRFMDTLGRRQPPYYITELFRAADKNKD
 NQICFDEFlyILGKLVKDYHLQFHRQLCAHYCTEHSly

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI



[View online »](#)

Cloning Scheme:


ACCN: NM_199422

ORF Size: 327 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.


RefSeq: [NM_199422.1](#), [NP_955454.1](#)

RefSeq Size: 741 bp

RefSeq ORF: 327 bp

Locus ID: 381493
 UniProt ID: [Q6S5I3](#)
 Cytogenetics: 3 F1
 MW: 12.9 kDa

Product images:



Circular map for MR220254