

Product datasheet for **MR220260L3V**

Spink6 (NM_001013797) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Spink6 (NM_001013797) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Spink6
Synonyms:	EG433180
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001013797
ORF Size:	315 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR220260).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001013797.1
RefSeq Size:	587 bp
RefSeq ORF:	318 bp
Locus ID:	433180
UniProt ID:	Q8BT20
Cytogenetics:	18 B3
Gene Summary:	Serine protease inhibitor selective for kallikreins. Efficiently inhibits KLK5 and human KLK2, KLK4, KLK5, KLK6, KLK7, KLK12, KLK13 and KLK14. Doesn't inhibit human KLK1 and KLK8. [UniProtKB/Swiss-Prot Function]



[View online »](#)