

Product datasheet for **MR220302**

Ccdc27 (NM_001033455) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Ccdc27 (NM_001033455) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Ccdc27
Synonyms:	Gm1035
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide
Sequence:

>MR220302 representing NM_001033455
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGAAGCAGGAGAACTTGACTCAGGAATGCTGCTCCACATCCCCAAGATCAGCAAGGGCCTGATGGTGC
TTCAGAGCATAGCCAGCCGGGGCTGCGACACACAGGACTCAGAGAAGAAACCGCAAAGTACCCAGCAAAG
CCTCAAAAAGTCAAGCCAGGCGGCTGCCGATACTTCTATGACAAGGAGGACCAGATCAGGAACTCAGC
AAACACAACGGCTTTTTGTCTGAAATGGAAGATATGAGAAAAGCCTTCTCATGCGGCCTGGGTGCCCC
AGTTTAGTACCAGGACAACATCCATGTCTCATGTGGTCTGCCATTATGGTTGATCTTCCAAGGACCTG
CTCTGGAGTTTGAAGTTGACTGAGGACCATCCCTAGGCAGGCTGGGCTCTGCCAGCAGTGTGGATGGG
AGGGTCTTTCCATTACAGCAAGAGTGCTTGTGAGCTCAACTACCCACGGAAGAGGAGTGAGCCTTCAGACC
CGAGTCCCACCGCAGCCCCACTGTGGTCAAGAAGTCCCAGAGGACCAGAACGCCCTGGTATATCTCAGT
CATCCATGAGAAGGATCACTCTCTGCTCCTGATGGGAGAAGAGTCCAGCGTTTCTCCGAGATGGAGAGC
CAGATGCAAGAAGAAAGATCAGGAGATCCTGACACTCCAGAAGGAGAAGGAGGCCCTGAAGAAGCAGCTGA
AGAACCTGCTCAGAGGCAAGGCCACAGAACTCTTCAGCTTCCATCAAGATGGACCGGTCTTTGAGAC
GCCCTGAAGCTGGGGAGGATGAGTGTGCTGAAGACCATCTACAAGAAGAGGATGAGCTACAGCACTGG
ATGCAGATGCAGGAGGAATATAGCATGGCGGAGAGTAGCAAGGAGCTGCACGTGGAACAGGAAAGCGCCA
TAGAGGAAAAGAGCTCTGAGGGGCCACCCGAGGAGGCTGCAGCTGCCAACTGAGCAGGCCCAGCCAGAG
CAAGACTGAGACGCTGCTGGAGTGGGCCCCGAGGAGGAGGAGGAGGAGGAAGAGGAGGAGGTGGAGGGT
GATGAAGCAAAAGGAACCTGAGGAGGGAGAGATATTGGTGAATGAGGAAGAGGCATCGTGGGAGCTTCGGG
AGGACGAGGAATGTCACCCCAAGAGGCTCTACTCCATGACGGAGTCTTTGAGGAGGAGCTGATGGCTCA
GCTGGAGGAGTATGAACGGATGTTGATGGACTTCCAGAGAGAGCTGGAATTCAGTCCAGATATTCC
CTGGCCACAGGGACCATCACATCTTTCAGCGGCAAACTGACTTCCAAGAGTCTCAGCTGCGGAAGGTCA
CCACGGAGAACGAGCTGCTGGAGAAGGAGCTCCGGGAACGGAAGCAGCAGATACAAGACATGACAGACAA
GTTCTCCAACCTCCGGGAGGAAAAGAAGCATCAGGAGATAATGGGACTCATTGAGAAGGAAAACCTGGTC
CTCCGACAGCAAGTTGCAGACTTGAAGATGGATCTCATCAGCAGTGAAGGACTATTAAGAGCTGAACA
CCCAAACCAAGAGCTGGAAGATCAGGTTAACACCGACAAGGACCACCTGAGGCGGTGGAAGGATTTGCA
TGACGACCTACAGACACGTAACGAGATCATCAACAGACAGGCAACAGACTCGTGTGGTCTGGAGGCC
ACCCAGGCCAGGTATGAGAAGCTAAGGAATAAGATCATCCAGGCTGTCTTCAGTGTGAGTGGGAACAAGA
ACTTGTCCATGGAACCTCTGATAGCTACATCCTGGAATCCCTGCAGAGGATCATCTCGGAGAGGAGTGA
CTTCTACAGTCAGCTGAAGCAAAAGGGCGTGAAGGTGCCCCACTACAGCAGTCAGATGTCTCCCTGCCT
AGCAAGATCAAGAAGATGGCATCCAAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR220302 representing NM_001033455
Red=Cloning site Green=Tags(s)


MKQENLTQECSTSPKISKGLMVLQSIASRGCDDTQDSEKKPQSTQQSLKKSSQAAAGYFYDKEDQIRKLS
 KHNGFLSEMEDMRKAFLMRPGCPQFSTRITSMHVGSAIMVDLPRTC SGVWKL TEDHPLGRLGSASSVDG
 RVFPFSK SACLNYPRKRSESPDPSPTGSP TVVKKSQRTRTPWYI SVIHEKDHSLLLMGEELQRFSEMS
 QMQKKDQEILTLQKEKALKQLKNLLRGKGTETSSASIKMDRSFETPLKLRMSVLKTIYKEEDELQHW
 MQMQEEYSMAESSKELHVEPGSAIEEKSSGEPPEEAAAALSRPSQSKTETLLEVPEEEEEEEEEVEG
 DEAKGTEEGEILVNEEEASWELREDEECHPKRSYSMTESFEEELMAQLEEYERMLMDFQRELEFTRSRYS
 LATGTITSLQRQDFQESQLRKVTTENELLEKELRERKQQIQDMDKFSNLREEKKHQEIMGLIEKENLV
 LRQQVADLKMDL ISSERTIKELNTQKLELDQVNTDKDHLRRWKDLHDDLQTRNEIIQQTEQQTRVVLEA
 TQARYEKLRNKIIQAVFSVSGNKNLSMELSDSYILESLQRIISERSDFYSQLKQKGVKVPPLQQSDVSLP
 SKIKKMASK

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN: NM_001033455

ORF Size: 1917 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001033455.2](#), [NP_001028627.1](#)

RefSeq Size: 2084 bp

RefSeq ORF: 1920 bp


Locus ID: 381580

UniProt ID: [Q3V036](#)

Cytogenetics: 4 E2

MW: 73.7 kDa

Product images:



Circular map for MR220302