

Product datasheet for **MR220308L4V**

Asb9 (NM_027027) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Asb9 (NM_027027) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Asb9
Synonyms:	1700011M07Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_027027
ORF Size:	873 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR220308).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_027027.2 , NP_081303.1
RefSeq Size:	1128 bp
RefSeq ORF:	873 bp
Locus ID:	69299
UniProt ID:	Q91ZT8
Cytogenetics:	X F5


[View online »](#)

Gene Summary:

Substrate-recognition component of a SCF-like ECS (Elongin-Cullin-SOCS-box protein) E3 ubiquitin-protein ligase complex which mediates the ubiquitination and subsequent proteasomal degradation of target proteins. Recognizes at least two forms of creatine kinase, CKB and CKMT1A (By similarity).[UniProtKB/Swiss-Prot Function]