

Product datasheet for MR220315

Hnrnpc (NM_016884) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Hnrnpc (NM_016884) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Hnrnpc
Synonyms:	AL022939; D14Wsu171e; hnrnp-C; hnRNP1; hnRNP2; Hnrpc; Hnrpc1; Hnrpc2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>MR220315 representing NM_016884 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCTAGCAATGTTACCAACAAGACAGATCCTCGGTCCATGAATCCCCTGTATTCATTGGGAATCTCA
ACACTCTGGTGGTCAAGAAGTCTGATGTGGAGGCCATCTTTCAAAGTATGGCAAATTTGGGGCTGCTC
TGTGCATAAAGGCTTTGCCTTTGTCCAGTATGTTAATGAAAGAAATGCCCGAGCTGCTGTAGCTGGCGAG
GATGGCAGAATGATTGCTGGCCAGGTTTTAGATATTAACCTGGCTGCAGAGCCAAAAGTGAACCGAGGAA
AAGCAGGTGTGAAACGATCTGCAGCGGAGATGTACGGTTCAGTACCAGAACACCCTTCTCCGTCCCCTCT
ACTCAGTTCCTCATTGACTTGGACTATGACTTTCAACGGGATTATTATGACAGGATGTACAGTTACCCA
GCGCGGGTTCTCCTCCTCCTCCTATTGCTCGAGCTGTGGTGCCTTCAAACGTCAGCGTGTTCAGGGA
ACACCTCACGAAGGGGCAAAAGTGGATTCAATTCGAAGAGTGGACAAAGGGGATCTTCTCCAAGTCTGG
AAAATTGAAAGGTGATGACCTTCAGGCCATTAAGGAGCTGACTCAGATAAAACAAAAGTGGATTCT
CTTCTGAAAAGCCTGGAAAAATGAAAAAGAACAAGCAAGCAAGCAGACTTGTCTTCTCATCCCCAG
TAGAGATGAAGAATGAAAAGTCTGAAGAAGAGCAGAGCAGCGCCTCTGTGAAGAAAGATGAGACTAATGT
GAAGATGGAGTCTGAGGCAGGTGCAGATGACTCTGCTGAGGAGGGTGACCTGCTGGATGATGACGATAAT
GAAGATCGGGGGATGACCAGCTGGAGTTGAAGGATGATGAAAAAGAGCCTGAGGAAGGAGAAGACGACA
GAGACAGCGCCAATGGGGAGGATGACTCT

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >MR220315 representing NM_016884
Red=Cloning site Green=Tags(s)

MASNVTNKTDPRSMNSRVFIGNLNLTLVVKKSDVEAIFSKYGKIVGCSVHKGFVQYVNERNARAAVAGE
 DGRMIAGQVLDINLAAEPKVNVRGKAGVKRSAAEMYGSVPEHPSPSPLSSSFDLDYDFQRDYYDRMYSYP
 ARVPPPPPIARAVVPSKRQRVSGNTSRRGKSGFNSKSGQRGSSSKSGKLGDDLQAIIKELTQIKQKVD
 LLESLKIEKEQSKQADLSFSSPVEMKNEKSEEEQSSASVKKDETNVKMESEAGADDSAEEDLLDDDDN
 EDRGDDQLLELKDDEKEPEEGEDDRDSANGEDDS

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/ja1225_e02.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF

ACCN: NM_016884

ORF Size: 939 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_016884.4](#)

RefSeq Size: 2844 bp

RefSeq ORF: 942 bp

Locus ID: 15381


UniProt ID: [Q9Z204](#)

Cytogenetics: 14 26.79 cM

MW: 34.8 kDa

Gene Summary: Binds pre-mRNA and nucleates the assembly of 40S hnRNP particles. Interacts with poly-U tracts in the 3' UTR or 5'-UTR of mRNA and modulates the stability and the level of translation of bound mRNA molecules. Single HNRNPC tetramers bind 230-240 nucleotides. Trimers of HNRNPC tetramers bind 700 nucleotides. May play a role in the early steps of spliceosome assembly and pre-mRNA splicing. N6-methyladenosine (m6A) has been shown to alter the local structure in mRNAs and long non-coding RNAs (lncRNAs) via a mechanism named 'm(6)A-switch', facilitating binding of HNRNPC, leading to regulation of mRNA splicing. [UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR220315