

Product datasheet for **MR220344L3V**

Defb41 (NM_183124) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Defb41 (NM_183124) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Defb41
Synonyms:	9230102D03Rik; BD-17; Bd-41; Defb16; Defb17; Gm15386; OTTMUSG00000020946; OTTMUSG00000020948
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_183124
ORF Size:	186 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR220344).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_183124.3 , NP_898947.1
RefSeq Size:	467 bp
RefSeq ORF:	189 bp
Locus ID:	77673
UniProt ID:	Q30KP6
Cytogenetics:	1 A3
Gene Summary:	Has bactericidal activity.[UniProtKB/Swiss-Prot Function]



[View online »](#)