

Product datasheet for **MR220351L3V**

Dpysl3 (NM_001136086) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Dpysl3 (NM_001136086) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Dpysl3
Synonyms:	CRMP; CRMP-4; DRP-3; TUC4; U; Ul; Ulip; ULIP-1; Ulip1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001136086
ORF Size:	1704 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR220351).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001136086.2 , NP_001129558.1
RefSeq Size:	5252 bp
RefSeq ORF:	1707 bp
Locus ID:	22240
Cytogenetics:	18 B3



[View online »](#)

Gene Summary:

This gene encodes a protein that belongs to the TUC (TOAD-64/Ulip/CRMP) family of proteins. Members of this family are phosphoproteins that function in axonal guidance and neuronal differentiation during development and regeneration of the nervous system. A mutation in the human gene is associated with amyotrophic lateral sclerosis. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Apr 2014]