

Product datasheet for **MR220353**

II20 (NM_021380) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: II20 (NM_021380) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: II20
Synonyms: Zcyto10
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >MR220353 representing NM_021380
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAAAGGCTTTGGTCTTGCCTTTGGACTGTTCTCCGCTGTGGGTTTTCTTCTCTGGACTCCTTTAACTG
GGCTCAAGACCCATTTGGGAAGCTGTGTGATTACTGCAAACCTACAGGCAATACAAAAGGAATTTTC
TGAGATTCGGGATAGTGTGCAAGCTGAAGATACAAATATTGACATCAGAATTTAAAGGACGACTGAGTCT
TTGAAAGACATAAAGTCTTTGGATAGGTGCTGCTTCTTCGTCATCTAGTGAGATTCTATCTGGACAGGG
TATTCAAAGTCTACCAGACCCCTGACCACCATACCCTGAGAAAGATCAGCAGCCTCGCCAACTCCTTTCT
TATCATCAAGAAGGACCTCTCAGTCTGTCACTTCTCACATGGCATGTCAATTGTGGGGAAGAAGCAATGGAG
AAATACAACCAAATTTCTGAGTCACTTCATAGAGTTGGAACCTTCAGGCAGCGGTGTTAAAGGCTTTGGGAG
AACTAGGCATTCTTCTGAGATGGATGGAGGAGATGCTA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR220353 representing NM_021380
Red=Cloning site Green=Tags(s)

MKGFGLAFGLFSAVGFLWLTPLTGLKTLHLGSCVITANLQAIQKEFSEIRDSVQAEDTNIDIRILRTTES
LKDIKSLDRCCFLRHLVRFYLDRVFKVYQTPDHHTLRKISSLANSLI IKKDL SVCHSHMACHGEEAME
KYNQILSHFIELELQAAVVKALGELGILLRWMEEML

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mm9006_h04.zip



[View online »](#)

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_021380

ORF Size: 528 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_021380.2](#), [NP_067355.1](#)

RefSeq Size: 824 bp

RefSeq ORF: 531 bp

Locus ID: 58181


UniProt ID: [Q9JKV9](#)

Cytogenetics: 1 E4

MW: 20.5 kDa

Gene Summary: Cytokine that may be involved in epidermal function and psoriasis. Acts through STAT3. [UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR220353