

Product datasheet for MR220382L3

Gstm3 (NM_010359) Mouse Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Gstm3 (NM_010359) Mouse Tagged Lenti ORF Clone
Tag:	Myc-DDK
Symbol:	Gstm3
Synonyms:	Fsc2; mGSTM5
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR220382).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF.

ACCN:	NM_010359
ORF Size:	654 bp



[View online »](#)

