

Product datasheet for **MR220402L4V**

Ctsa (NM_008906) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Ctsa (NM_008906) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Ctsa
Synonyms:	AU019505; P; Pp; PPCA; Ppqb
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_008906
ORF Size:	1476 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR220402).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_008906.3 , NP_032932.2
RefSeq Size:	3674 bp
RefSeq ORF:	1479 bp
Locus ID:	19025
Cytogenetics:	2 85.27 cM



[View online »](#)

Gene Summary:

This gene encodes a glycoprotein with deamidase, esterase and carboxypeptidase activities. The encoded protein associates with and provides a protective function to the lysosomal enzymes beta-galactosidase and neuraminidase. Deficiency of the related gene in humans results in galactosialidosis. The proprotein is processed into two shorter chains. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Jan 2013]