

## Product datasheet for **MR220406L3V**

### Katnal2 (NM\_027721) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Katnal2 (NM_027721) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Katnal2
Synonyms:	3110023G01Rik; 4933439B08Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_027721
ORF Size:	1230 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR220406).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_027721.2</a> , <a href="#">NP_081997.1</a>
RefSeq Size:	1673 bp
RefSeq ORF:	1230 bp
Locus ID:	71206
UniProt ID:	<a href="#">Q9D3R6</a>
Cytogenetics:	18 E3
Gene Summary:	Severs microtubules in vitro in an ATP-dependent manner. This activity may promote rapid reorganization of cellular microtubule arrays.[UniProtKB/Swiss-Prot Function]



[View online »](#)