

## Product datasheet for MR220431L4

### Bmp2k (NM\_080708) Mouse Tagged Lenti ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Bmp2k (NM_080708) Mouse Tagged Lenti ORF Clone
Tag:	mGFP
Symbol:	Bmp2k
Synonyms:	4933417M22Rik; AA673486; AV128808; BIKE
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR220431).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



\* The last codon before the Stop codon of the ORF.

ACCN:	NM_080708
ORF Size:	3414 bp



[View online »](#)

