

Product datasheet for MR220456L4

Adamtsl4 (NM_144899) Mouse Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Adamtsl4 (NM_144899) Mouse Tagged Lenti ORF Clone
Tag:	mGFP
Symbol:	Adamtsl4
Synonyms:	ADAMTSL-4; Tsr; Tsrc1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR220456).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF.


ACCN:	NM_144899
ORF Size:	3111 bp



[View online »](#)

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_144899.2 , NP_659148.2
RefSeq Size:	3938 bp
RefSeq ORF:	3111 bp
Locus ID:	229595
UniProt ID:	Q80T21
Cytogenetics:	3 F2.1
Gene Summary:	The protein encoded by this gene is a member of the ADAMTS superfamily of secreted proteins, which contain a metalloprotease domain at the N-terminus and a C-terminal ancillary domain. ADAMTS-like proteins lack protease activity and resemble the ancillary domain of ADAMTS proteins. ADAMTS-like proteins have been implicated in regulation of the extracellular matrix. The encoded protein contains 7 thrombospondin type 1 repeats, a conserved extracellular domain. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Sep 2014]

Product images:



Circular map for MR220456L4