

Product datasheet for **MR220459L4V**

Erp27 (NM_026983) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Erp27 (NM_026983) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Erp27
Synonyms:	1810033M07Rik; 1810047B09Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_026983
ORF Size:	816 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR220459).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_026983.2 , NP_081259.1
RefSeq Size:	953 bp
RefSeq ORF:	819 bp
Locus ID:	69187
UniProt ID:	Q9D8U3
Cytogenetics:	6 G1
Gene Summary:	Specifically binds unfolded proteins and may recruit protein disulfide isomerase PDIA3 to unfolded substrates. Binds protein substrates via a hydrophobic pocket in the C-terminal domain. May play a role in the unfolded stress response.[UniProtKB/Swiss-Prot Function]



[View online »](#)