

Product datasheet for **MR220480L4V**

Etfb (NM_026695) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Etfb (NM_026695) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Etfb
Synonyms:	0610009I16Rik; 2810441H06Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_026695
ORF Size:	768 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR220480).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_026695.2
RefSeq Size:	920 bp
RefSeq ORF:	768 bp
Locus ID:	110826
UniProt ID:	Q9DCW4
Cytogenetics:	7 28.25 cM



[View online »](#)

Gene Summary:

Heterodimeric electron transfer flavoprotein that accepts electrons from several mitochondrial dehydrogenases, including acyl-CoA dehydrogenases, glutaryl-CoA and sarcosine dehydrogenase. It transfers the electrons to the main mitochondrial respiratory chain via ETF-ubiquinone oxidoreductase (By similarity). Required for normal mitochondrial fatty acid oxidation and normal amino acid metabolism (PubMed:25023281). ETFB binds an AMP molecule that probably has a purely structural role (By similarity).[UniProtKB/Swiss-Prot Function]